

Standard	Dilution	Contact time	Temperature	Strains
EN 1040 - Bactericidal	ready to use	15 minutes	20°C	<i>Pseudomonas aeruginosa, Staphylococcus aureus</i>
EN 1276 - Bactericidal	ready to use	5 minutes	20°C	<i>Pseudomonas aeruginosa, Staphylococcus aureus, Enterococcus hirae, Escherichia Coli</i>
EN 1276 - Additional Bactericidal	ready to use	5 minutes	20°C	<i>Salmonella enterica subsp. salamae serotyoe Tranoroa, Listeria monocytogene, Legionella pneumophila</i>
EN 13697 - Bactericidal	ready to use	5 minutes	20°C	<i>Pseudomonas aeruginosa, Staphylococcus aureus, Enterococcus hirae, Escherichia Coli</i>
EN 13697 - Fungicidal	ready to use	15 minutes	20°C	<i>Candida Albicans, Aspergillus niger</i>

Material to be used exclusively by the Sutter Professional sales force